



PHARMACEUTICAL EQUIPMENT AND SERVICES, L.L.C.



- Removes cups and punch tips from unwanted tableting punches safely and efficiently
- Allows hands-free destruction of pharmaceutical tooling without generating heat, dust, or noise
- Constructed completely from stainless steel
- Designed for installation in existing tooling rooms or any typical pharmaceutical environment

**ERASER 2620® ...
Safe, Automated
Destruction of
Unwanted Tablet
Press Tooling**



UNIT LOADED WITH B TOOLING

**ERASER 2620®
PATENTED**

ERASER 2620[®]
PATENTED

Dear Customer,

Safe, simple destruction of unwanted punches is required by all users of compression tooling. Until now, most companies destroy tooling by grinding, hammering, or torching the punch tips. These practices are dangerous and require some level of skill to perform without injury. I have developed the ERASER 2620[®] as a user-friendly device that safely and effectively destroys unwanted punches. This self-contained unit installed in the tooling operations area is ideal at its task, and gives closure to issues involving the care of unwanted punches. I am willing to send a unit or invite you to visit our facility in New Jersey for a demonstration. These units come with a 2-year warranty for parts and service. Please contact us at your convenience, and thank you for your support and consideration.

Sincerely,



Steven McLean
General Manager

DIMENSIONS

Width	33"
Depth	43"
Height	58"
Weight	402 lbs
Voltage Requirements	110 VAC 20 amps



PHARMACEUTICAL EQUIPMENT AND SERVICES, L.L.C.

5042 Industrial Road, Unit D • Wall, NJ 07719
732-714-2250 • fax 732-714-2255
www.stevens-industries.com